

Laboratory Mixing Molder

LMM



Description

The versatile Laboratory Mixing Molder is designed to mold miniature specimens for dynamic, tensile and impact testing. This benchtop mixer, extruder and molding machine will produce cost-effective moldings using less than one gram of material.

This compact R & D tool allows production of cost-effective moldings for testing standard or newly developed thermoplastic material.

The LMM prepares a polymer melt for injection using two types of mixing: extensive and intensive. This combination ensures material homogeneity for the feed process.

Features

- Molds miniature specimen for dynamic, tensile and impact
- Effective moldings from less than one gram of material
- Interchangeable cup capacities
- Sample cycles in less than 30 seconds
- Digital display PID temperature controller
- Operating temperature from ambient to 400°C
- Standard C-clamp to secure molds
- Accommodates mold lengths of 1" to 4" (2.54cm to 10cm)
- Variable speed rotor control

Optional

- Heated C-clamp, with integral temperature controller to optimize mold temperature to process a wide variety of geometrical shapes

Intensive mixing produces high shear rates to break up clumps of material. Extensive mixing causes a folding type of material action to uniformly distribute the various ingredients of the melt. Mixing begins when the sample materials are heated through the conduction and mechanical shearing of the rotor. From mixing to a finished mold, the LMM is capable of cycle times less than 30 seconds.

www.dynisco.com

Dynisco
38 Forge Parkway
Franklin, MA 02038
USA

Hotline 1-800-DYNISCO
Phone +1-508-541-9400
Fax +1-508-541-6206
Email infoinst@dynisco.com

Dynisco Europe GmbH
Pfaffenstr. 21
74078 Heilbronn
Germany

Phone +49 7131 297-0
Fax +49 7131 23260
Email dyniscoeurope@dynisco.com

Procedure

Once the processing temperature is stabilized, the sample materials are loaded into the stator cup area. The rotor is then lowered into the cup and mixing begins. Once the material is melted, the radial mixing step begins as the rotor

is cycled up and down. The pressure from both the rotor and the elastic melt extruder passes the fully melted and mixed polymer into the specimen mold. The process is completed with the removal of the specimen from the mold.

Specifications

Dimensions	27"W x 12"D x 24"H (68cm x 30cm x 61cm)
Standard Specimen Size	Cylindrical dumbbell Length: .875" (22.2mm) Diameter: .06" (1.6mm) Weight: .02g
Cup Capacities	2.0cc or 4.0cc
Injection Force	0.5kg force applied to injection lever exerts a downward force of 5.3kg or 4.2kg/cm ²
Operating Temperature Range	Ambient to 450°C
Weight	(30kg) 65 lbs.
Electrical	120/240V, 50/60Hz, 4A

Specimen Types

- Bar
- Coil
- Cylindrical dumbbell
- Tensile and impact
- Surgical implants
- User specified

Stocked, Distributed, and Supported by

SENSORS
INCORPORATED

507 Kelsey Street • Delano, MN 55328
Phone 763-972-1040 Fax 763-972-1041
Toll Free 888-920-0939
Sensorsincorporated.com

All dimensions are in inches (millimeters).

©2008. Dynisco reserves the right to make changes without notice.
Refer to www.dynisco.com for access to Operator Manual and other support documentation.
DDS276070 10.08